

LHB Polyclonal Antibody



Catalog Number: 163972

Product name

LHB Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to LHB

Preparation

Antigen: A synthetic peptide of human LHB

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100

Validations



Western blot - LHB Polyclonal Antibody

Western blot analysis of extracts of K-562 cells, using LHB antibody at 1:400 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.