

MRP1 / ABCC1 Polyclonal Antibody



Catalog Number: 163985

Product name

MRP1 / ABCC1 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to MRP1 / ABCC1

Preparation

Antigen: A synthetic peptide of human MRP1 / ABCC1

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200

Validations



Western blot - MRP1 / ABCC1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using MRP1 / ABCC1 antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.