

MMP10 Polyclonal Antibody



Catalog Number: 163986

Product name

MMP10 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to MMP10

Preparation

Antigen: A synthetic peptide of human MMP10

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations




Western blot - MMP10 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using MMP10 antibody at 1:400 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunohistochemistry - MMP10 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat kidney using MMP10 antibody at dilution of 1:100



MMP10 Polyclonal Antibody



Catalog Number: 163986

Immunohistochemistry - MMP10 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using MMP10 antibody at dilution of 1:100 (40x lens).