

# TCF7 Polyclonal Antibody



Catalog Number: 164015

## Product name

TCF7 Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to TCF7

## Preparation

Antigen: A synthetic peptide of human TCF7

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:1000

## Validations



Western blot - TCF7 Polyclonal Antibody

Western blot analysis of extracts of mouse spleen, using TCF7 antibody at 1:300 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.