

ADCYAP1R1 Polyclonal Antibody



Catalog Number: 164029

Product name

Rabbit IgG

ADCYAP1R1 Polyclonal Antibody

Applications

Specificity

WB

Human

Dilutions

Antibody description


WB 1:500 - 1:2000

Polyclonal antibody to ADCYAP1R1

Validations

Preparation

Antigen: A synthetic peptide of human ADCYAP1R1



Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Western blot - ADCYAP1R1 Polyclonal Antibody

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Western blot analysis of extracts of various cell lines, using ADCYAP1R1 antibody at 1:400 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.

Clonality

Polyclonal

Ig Type