

PPAR alpha Polyclonal Antibody



Catalog Number: 164032

Product name

Polyclonal

PPAR alpha Polyclonal Antibody

Ig Type**Specificity**

Rabbit IgG

Human, Mouse, Rat

Applications**Antibody description**

WB

Polyclonal antibody to PPAR alpha

Dilutions**Preparation**

WB 1:500 - 1:1000

Antigen: A synthetic peptide of human PPAR alpha

Validations**Formulation**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Western blot - PPAR alpha Polyclonal Antibody

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Western blot analysis of extracts of 231 cells, using PPAR alpha antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.

Clonality