

SSTR5 Polyclonal Antibody



Catalog Number: 164037

Product name

SSTR5 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to SSTR5

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 275-364 of human SSTR5 (NP_001044.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:1000

Validations



Western blot - SSTR5 Polyclonal Antibody

Western blot analysis of extracts of HT-29 cells, using SSTR5 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.