

SLC22A6 Polyclonal Antibody



Catalog Number: 164054

Product name

SLC22A6 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to SLC22A6

Preparation

Antigen: A synthetic peptide of human SLC22A6

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - SLC22A6 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using SLC22A6 antibody at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.