

ABCA6 Polyclonal Antibody



Catalog Number: 164056

Product name

ABCA6 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to ABCA6

Preparation

Antigen: Recombinant protein of mouse ABCA6

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations



Western blot - ABCA6 Polyclonal Antibody

Western blot analysis of extracts of 293T cells, using ABCA6 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunohistochemistry - ABCA6 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse liver using ABCA6 antibody at dilution of 1:100 (40x lens).