

FAM160B2 Polyclonal Antibody



Catalog Number: 164057

Product name

FAM160B2 Polyclonal Antibody

WB 1:500 - 1:2000

IHC 1:50 - 1:200

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to FAM160B2

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 420-520 of human FAM160B2 (NP_073586.5).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

WB, IHC


Dilutions

Validations



Immunohistochemistry - FAM160B2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse liver fibrosis using FAM160B2 antibody at dilution of 1:200 (40x lens).



Immunohistochemistry - FAM160B2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse heart using FAM160B2 antibody at dilution of 1:200 (40x lens).