

iNOS Polyclonal Antibody



Catalog Number: 164061

Product name

iNOS Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to iNOS

Preparation

Antigen: Recombinant protein of human iNOS

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF, FC

Dilutions


WB 1:500 - 1:1000

IHC 1:50 - 1:100

IF 1:20 - 1:50

FC 1:20 - 1:50

Validations



Western blot - iNOS Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using iNOS antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.