

SOX1 Polyclonal Antibody



Catalog Number: 164068

Product name

SOX1 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to SOX1

Preparation

Antigen: A synthetic peptide of human SOX1

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, ICC


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200

ICC 1:50 - 1:200

Validations



Western blot - SOX1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using SOX1 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.