

# COX8A Polyclonal Antibody



Catalog Number: 164081

## Product name

COX8A Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to COX8A

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-69 of human COX8A (NP\_004065.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:1000 - 1:2000

## Validations



Western blot - COX8A Polyclonal Antibody

Western blot analysis of extracts of mouse skeletal muscle, using COX8A antibody .