

NFIB Polyclonal Antibody



Catalog Number: 164088

Product name

NFIB Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to NFIB

Preparation

Antigen: Recombinant protein of human NFIB

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:1000

Validations



Western blot - NFIB Polyclonal Antibody

Western blot analysis of extracts of human heart, using NFIB antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.