

SEC23IP Polyclonal Antibody



Catalog Number: 164106

Product name

SEC23IP Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to SEC23IP

Preparation

Antigen: A synthetic peptide of human SEC23IP

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - SEC23IP Polyclonal Antibody

Western blot analysis of extracts of HeLa cells, using SEC23IP antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.