

FAM120A Monoclonal Antibody



Catalog Number: 164108

Product name

FAM120A Monoclonal Antibody

Specificity

Human

Antibody description

Monoclonal antibody to FAM120A

Preparation

Antigen: A synthetic peptide of human FAM120A

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:1000

IHC 1:20 - 1:50

IF 1:20 - 1:50

Validations



Western blot - FAM120A Monoclonal Antibody

Western blot analysis of extracts of various cell lines, using FAM120A antibody. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.