

# ZNF703 Polyclonal Antibody



Catalog Number: 164134

## Product name

ZNF703 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to ZNF703

## Preparation

Antigen: A synthetic peptide of human ZNF703

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:1000 - 1:2000

## Validations



Western blot - ZNF703 Polyclonal Antibody

Western blot analysis of extracts of Jurkat cells, using ZNF703 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.