

SLC25A25 Polyclonal Antibody



Catalog Number: 164139

Product name

SLC25A25 Polyclonal Antibody

Specificity

Human, Rat

Antibody description

Polyclonal antibody to SLC25A25

Preparation

Antigen: A synthetic peptide of human SLC25A25

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:1000

IHC 1:50 - 1:100

Validations



Western blot - SLC25A25 Polyclonal Antibody

Western blot analysis of extracts of U-937 cells, using SLC25A25 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.