

NOSTRIN Polyclonal Antibody



Catalog Number: 164140

Product name

NOSTRIN Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to NOSTRIN

Preparation

Antigen: A synthetic peptide of human NOSTRIN

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:1000

Validations



Western blot - NOSTRIN Polyclonal Antibody

Western blot analysis of extracts of HepG2 cells, using NOSTRIN antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.