

DOCK2 Polyclonal Antibody



Catalog Number: 164151

Product name

DOCK2 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to DOCK2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1551-1830 of human DOCK2 (NP_004937.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:1000

Validations



Western blot - DOCK2 Polyclonal Antibody

Western blot analysis of extracts of Jurkat cells, using DOCK2 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.