

# DACH1 Polyclonal Antibody



Catalog Number: 164211

---

## Product name

DACH1 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to DACH1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 450-570 of human DACH1 (NP\_542937.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:200 - 1:2000

IHC 1:20 - 1:200

## Validations



Immunohistochemistry - DACH1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using DACH1 antibody at dilution of 1:100 (40x lens).