

DUSP4 Polyclonal Antibody



Catalog Number: 164218

Product name

DUSP4 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to DUSP4

Preparation

Antigen: Recombinant protein of human DUSP4

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:200 - 1:2000

IHC 1:50 - 1:200

Validations



Western blot - DUSP4 Polyclonal Antibody

Western blot analysis of extracts of A-549 cells, using DUSP4 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.