

FUT6 Polyclonal Antibody

Catalog Number: 164225

Product name

FUT6 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to FUT6

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 170-359 of human FUT6 (NP_000141.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:100 - 1:200

Validations




Western blot - FUT6 Polyclonal Antibody

Western blot analysis of extracts of SW480 cells, using FUT6 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time:



Immunohistochemistry - FUT6 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver using FUT6 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - FUT6 Polyclonal Antibody

FUT6 Polyclonal Antibody



Catalog Number: 164225

Immunohistochemistry of paraffin-embedded

human stomach using FUT6 antibody at dilution of 1:100 (40x lens).