

SMAD9 Polyclonal Antibody



Catalog Number: 164244

Product name

SMAD9 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to SMAD9

Preparation

Antigen: A synthetic peptide of human SMAD9

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:1000

IF 1:20 - 1:50

Validations



Western blot - SMAD9 Polyclonal Antibody

Western blot analysis of extracts of A-549 cells, using SMAD9 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.