

# SEPT5 Polyclonal Antibody



Catalog Number: 164257

## Product name

SEPT5 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to SEPT5

## Preparation

Antigen: A synthetic Peptide of human SEPT5

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC


## Dilutions

WB 1:500 - 1:1000

IHC 1:20 - 1:50

## Validations

Human brain



Western blot - SEPT5 Polyclonal Antibody

Western blot analysis of extracts of human brain, using SEPT5 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.