

TRAPPC2 Polyclonal Antibody



Catalog Number: 164270

Product name

TRAPPC2 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to TRAPPC2

Preparation

Antigen: A synthetic Peptide of human TRAPPC2

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:1000

Validations



Western blot - TRAPPC2 Polyclonal Antibody

Western blot analysis of extracts of NCI-H292 cells, using TRAPPC2 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.