

EEF1E1 Monoclonal Antibody



Catalog Number: 164306

Product name

EEF1E1 Monoclonal Antibody

Specificity

Human, Mouse

Antibody description

Monoclonal antibody to EEF1E1

Preparation

Antigen: A synthetic peptide of human EEF1E1

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG

Applications

WB, IF, FC


Dilutions

WB 1:500 - 1:1000

IF 1:20 - 1:50

FC 1:20 - 1:50

Validations



Western blot - EEF1E1 Monoclonal Antibody

Western blot analysis of extracts of HepG2 cells, using EEF1E1 antibody. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.