

LMOD1 Polyclonal Antibody



Catalog Number: 164367

Product name

LMOD1 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to LMOD1

Preparation

Antigen: A synthetic Peptide of human LMOD1

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, FC


Dilutions

WB 1:500 - 1:1000

IHC 1:20 - 1:50

FC 1:20 - 1:50

Validations



Western blot - LMOD1 Polyclonal Antibody

Western blot analysis of extracts of HeLa cells, using LMOD1 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.