

# PHF11 Polyclonal Antibody



Catalog Number: 164382

## Product name

PHF11 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to PHF11

## Preparation

Antigen: A synthetic Peptide of human PHF11

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:1000

## Validations



Western blot - PHF11 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using PHF11 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.