

CSGALNACT2 Polyclonal Antibody



Catalog Number: 164403

Product name

CSGALNACT2 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to CSGALNACT2

Preparation

Antigen: A synthetic Peptide of human CSGALNACT2

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, FC


Dilutions

WB 1:500 - 1:1000

IHC 1:20 - 1:50

FC 1:20 - 1:50

Validations



Western blot - CSGALNACT2 Polyclonal Antibody

Western blot analysis of extracts of 293T cells, using CSGALNACT2 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.