

IWS1 Polyclonal Antibody



Catalog Number: 164405

Product name

IWS1 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to IWS1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 570-819 of human IWS1 (NP_060439.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:200

Validations



Immunofluorescence - IWS1 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using IWS1 antibody .