

# FBRS Polyclonal Antibody



Catalog Number: 164418

## Product name

FBRS Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to FBRS

## Preparation

Antigen: Recombinant protein of human FBRS

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:200 - 1:2000

IHC 1:50 - 1:200

IF 1:50 - 1:100

## Validations



## Western blot - FBRS Polyclonal Antibody

Western blot analysis of extracts of Jurkat cells, using FBRS antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.