

PRKCDBP Polyclonal Antibody



Catalog Number: 164444

Product name

PRKCDBP Polyclonal Antibody

Specificity

Human, Rat

Antibody description

Polyclonal antibody to PRKCDBP

Preparation

Antigen: A synthetic peptide of human PRKCDBP

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:200 - 1:1000

Validations



Western blot - PRKCDBP Polyclonal Antibody

Western blot analysis of extracts of rat kidney, using PRKCDBP antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.