

# DHX36 Polyclonal Antibody



Catalog Number: 164447

## Product name

DHX36 Polyclonal Antibody

## Specificity

Mouse, Rat

## Antibody description

Polyclonal antibody to DHX36

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human DHX36 (NP\_065916.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IP

## Dilutions

WB 1:200 - 1:2000

IP 1:50 - 1:100

## Validations



## Western blot - DHX36 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using DHX36 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.