

APT_X Polyclonal Antibody



Catalog Number: 164490

Product name

APT_X Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to APT_X

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 93-342 of human APT_X (NP_001182178.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF, IP


Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:100


IP 1:20 - 1:50

Validations




Western blot - APT_X Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using APT_X antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunofluorescence - APT_X Polyclonal Antibody

Immunofluorescence analysis of A549 cells using



APT_X antibody .

Immunoprecipitation - APT_X Polyclonal Antibody

APT_X Polyclonal Antibody



Catalog Number: 164490

Immunoprecipitation analysis of 150ug extracts of A549 cells using 3ug APT_X antibody . Western blot was performed from the immunoprecipitate using APT_X antibody at a dilution of 1:500.