

EIF3E Polyclonal Antibody



Catalog Number: 164518

Product name

EIF3E Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to EIF3E

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-445 of human EIF3E (NP_001559.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:100

Validations



Western blot - EIF3E Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using EIF3E antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunofluorescence - EIF3E Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using EIF3E antibody .