

# CYP4A11 Polyclonal Antibody



Catalog Number: 164523

---

## Product name

CYP4A11 Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to CYP4A11

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human CYP4A11 (NP\_000769.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100

## Validations



Immunohistochemistry - CYP4A11 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using CYP4A11 antibody at dilution of 1:100 (40x lens).