

TRIM22 Polyclonal Antibody



Catalog Number: 164540

Product name

TRIM22 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to TRIM22

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 215-411 of human TRIM22 (NP_006065.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:1000


IF 1:50 - 1:100

Validations



Western blot - TRIM22 Polyclonal Antibody

Western blot analysis of extracts of A-549 cells, using TRIM22 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunofluorescence - TRIM22 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using TRIM22 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.