

PSCA Polyclonal Antibody



Catalog Number: 164567

Product name

PSCA Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to PSCA

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 12-114 of human PSCA (NP_005663.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - PSCA Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using PSCA antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit .