

DOK1 Polyclonal Antibody

Catalog Number: 164595

Product name

DOK1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to DOK1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human DOK1 (NP_001372.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:100

Validations



Western blot - DOK1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using DOK1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.




Immunohistochemistry - DOK1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human breast cancer using DOK1 Antibody at dilution of 1:100 (40x lens).

DOK1 Polyclonal Antibody




Catalog Number: 164595



Immunohistochemistry - DOK1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat brain using DOK1 Antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - DOK1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human breast using DOK1 Antibody at dilution of



1:100 (40x lens).


Immunohistochemistry - DOK1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human gastric cancer using DOK1 Antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - DOK1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse heart using DOK1 Antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - DOK1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse pancreas using DOK1 Antibody at dilution of 1:100 (40x lens).