

# PCGF6 Polyclonal Antibody



Catalog Number: 164627

## Product name

PCGF6 Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to PCGF6

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 161-350 of human PCGF6 (NP\_001011663.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:100

## Validations



Western blot - PCGF6 Polyclonal Antibody

Western blot analysis of extracts of MCF7 cells, using PCGF6 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.



Immunofluorescence - PCGF6 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using PCGF6 antibody.