

# TOPBP1 Polyclonal Antibody



Catalog Number: 164636

## Product name

TOPBP1 Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to TOPBP1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1253-1522 of human TOPBP1 (NP\_008958.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:100 - 1:200


IF 1:50 - 1:200

## Validations



Western blot - TOPBP1 Polyclonal Antibody


Western blot analysis of extracts of 293T cells, using TOPBP1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure



time: 50s.

Immunohistochemistry - TOPBP1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse heart using TOPBP1 antibody at dilution of



1:100 (40x lens).

# TOPBP1 Polyclonal Antibody



Catalog Number: 164636

---

Immunofluorescence - TOPBP1 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using TOPBP1 antibody . Blue: DAPI for nuclear staining.