

# CYP2J2 Polyclonal Antibody



Catalog Number: 164642

## Product name

CYP2J2 Polyclonal Antibody

## Specificity

Mouse, Rat

## Antibody description

Polyclonal antibody to CYP2J2

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 200-300 of human CYP2J2 (NP\_000766.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



## Western blot - CYP2J2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using CYP2J2 Antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.