

PRMT2 Polyclonal Antibody



Catalog Number: 164657

Product name

PRMT2 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to PRMT2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human PRMT2 (NP_996845.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:100

Validations



Western blot - PRMT2 Polyclonal Antibody

Western blot analysis of extracts of mouse lung, using PRMT2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure



time: 90s.


Immunohistochemistry - PRMT2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver injury using PRMT2 antibody at dilution of 1:100 (40x lens).

PRMT2 Polyclonal Antibody



Catalog Number: 164657



Immunofluorescence - PRMT2 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using PRMT2 antibody .