

# SF3A2 Polyclonal Antibody



Catalog Number: 164699

## Product name

SF3A2 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to SF3A2

## Preparation

Antigen: A synthetic peptide of human SF3A2

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



Western blot - SF3A2 Polyclonal Antibody

Western blot analysis of extracts of HeLa cells, using SF3A2 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.