

# DICER1 Polyclonal Antibody



Catalog Number: 164713

## Product name

DICER1 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to DICER1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1200-1460 of human DICER1 (NP\_085124.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB, IHC, IF

## Dilutions


WB 1:500 - 1:1000  
IHC 1:50 - 1:200  
IF 1:100 - 1:200

## Validations



Western blot - DICER1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using DICER1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - DICER1 Polyclonal Antibody



Immunohistochemistry of paraffin-embedded human lung cancer using DICER1 antibody at dilution of 1:100 (40x lens).

# DICER1 Polyclonal Antibody



Catalog Number: 164713

---



Immunohistochemistry - DICER1 Polyclonal Antibody

Immunohistochemistry - DICER1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human gastric cancer using DICER1 antibody at dilution of 1:100 (40x lens).

Immunohistochemistry of paraffin-embedded mouse brain using DICER1 antibody at dilution of 1:100 (40x lens).