

# TNRC6A Polyclonal Antibody



Catalog Number: 164729

## Product name

TNRC6A Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to TNRC6A

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 950-1050 of human TNRC6A (NP\_055309.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:1000


IF 1:50 - 1:100

## Validations



Western blot - TNRC6A Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using TNRC6A antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL



Basic Kit .Exposure time: 90s.

Immunofluorescence - TNRC6A Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using TNRC6A antibody at dilution of 1:100. Blue: DAPI for nuclear staining.