

TRPM2 Polyclonal Antibody



Catalog Number: 164739

Product name

TRPM2 Polyclonal Antibody

Specificity

Mouse, Rat

Antibody description

Polyclonal antibody to TRPM2

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 900-1000 of human TRPM2 (NP_003298.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions


IF 1:50 - 1:100

Validations



Immunofluorescence - TRPM2 Polyclonal Antibody

Immunofluorescence analysis of rat brain using TRPM2 antibody at dilution of 1:100. Blue: DAPI



for nuclear staining.

Immunofluorescence - TRPM2 Polyclonal Antibody

Immunofluorescence analysis of mouse brain using TRPM2 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.