

CTCFL Polyclonal Antibody



Catalog Number: 164746

Product name

CTCFL Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to CTCFL

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human CTCFL (NP_001255970.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:20 - 1:100

Validations



Western blot - CTCFL Polyclonal Antibody
Western blot analysis of extracts of various cell lines, using CTCFL antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.

Immunofluorescence - CTCFL Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using CTCFL antibody. Blue: DAPI for nuclear staining.