

DAP Polyclonal Antibody



Catalog Number: 164747

Product name

DAP Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to DAP

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-102 of human DAP (NP_004385.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF

Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200

IF 1:20 - 1:100

Validations



Western blot - DAP Polyclonal Antibody
Western blot analysis of extracts of various cell lines, using DAP antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.

Immunohistochemistry - DAP Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat kidney using DAP antibody at dilution of 1:200 (40x lens).

Immunofluorescence - DAP Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using DAP antibody. Blue: DAPI for nuclear staining.